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Private Praticice in Oral Surgery and Implantology,
Clínica Delille, Coimbra, Portugal

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1. Indication profile (CASE 1)

Region	<input type="checkbox"/> aesthetic region	<input checked="" type="checkbox"/> non-aesthetic region
	<input checked="" type="checkbox"/> single tooth gap	<input type="checkbox"/> multiple tooth gaps
	<input checked="" type="checkbox"/> tooth 26	
Bone situation	<input checked="" type="checkbox"/> bone defect present	<input type="checkbox"/> no bone defect present
	<input checked="" type="checkbox"/> maxillary atrophy	
Soft tissue situation	<input type="checkbox"/> recession	<input type="checkbox"/> no recession
	<input type="checkbox"/> inflamed	<input type="checkbox"/> infected
	<input checked="" type="checkbox"/> thick biotype	<input type="checkbox"/> thin biotype
	<input checked="" type="checkbox"/> primary wound closure possible	<input type="checkbox"/> primary wound closure not possible
	<input type="checkbox"/> intact papillae	<input type="checkbox"/> impaired, missing papillae
Implantation	<input checked="" type="checkbox"/> adequately keratinised mucosa	<input type="checkbox"/> inadequate keratinised mucosa
	<input checked="" type="checkbox"/> simultaneously with bone augmentation (1 step)	<input type="checkbox"/> uneventful
	<input type="checkbox"/> successively to bone augmentation (2 steps)	

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Further indication sheets

> To receive these by mail free of charge, please contact: www.geistlich.com/indicationsheets.
> If you do not wish to collect indication sheets any more, please unsubscribe via your local distribution partner.

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1. Indication profile (CASE 2)

Region	<input checked="" type="checkbox"/> aesthetic region	<input checked="" type="checkbox"/> non-aesthetic region	
Bone situation	<input type="checkbox"/> single tooth gap	<input checked="" type="checkbox"/> multiple tooth gaps	
	<input type="checkbox"/> bone defect present	<input type="checkbox"/> no bone defect present	
Soft tissue situation	<input checked="" type="checkbox"/> maxillary atrophy	<input type="checkbox"/>	
	<input type="checkbox"/> recession	<input type="checkbox"/> no recession	
	<input type="checkbox"/> inflamed	<input type="checkbox"/> infected	
	<input checked="" type="checkbox"/> thick biotype	<input type="checkbox"/> thin biotype	
Implantation	<input checked="" type="checkbox"/> primary wound closure possible	<input type="checkbox"/> primary wound closure not possible	
	<input type="checkbox"/> intact papillae	<input checked="" type="checkbox"/> impaired, missing papillae	
	<input type="checkbox"/> adequately keratinised mucosa	<input type="checkbox"/> inadequate keratinised mucosa	<input type="checkbox"/> uneventful
	<input checked="" type="checkbox"/> simultaneously with bone augmentation (1 step)	<input type="checkbox"/>	<input type="checkbox"/>
	<input checked="" type="checkbox"/> successively to bone augmentation (2 steps)	<input type="checkbox"/>	<input type="checkbox"/>

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©Geistlich Pharma AG